

Bombas de infusión SMARTeZ Infusiones de corta duración



Tabla del efecto de una mayor o menor cantidad de relleno

Nº de modelo (REF)	481030	481090	481150	481040	481120	481100	481010	481050	481080	481020	481060
Caudal (mL/h)	50	50	50	100	100	100	200	250	250	500	175
Volumen de llenado /reposición (mL)	50	100	250	100	250	400	100	250	500	250	250
Volumen mínimo (ml)	50	75	150	50	200	275	50	150	360	250	250
Volumen máximo (ml)	60	110	275	120	300	550	100	300	550	300	300
Volumen residual (ml)	2	2	3	2	3	5	2	3	5	3	3

Dauer der Infusion

Minutos	Volumen de llenado/reposición (mL)										
00:15							50(*)				
00:30				50			100	150		250	
00:45				75				200		300	140
01:00	50			100				250			175
01:15	60			120				300			220
01:30		75							360		250
01:45									400		300
02:00		100			200				475		
02:15		110			220				500		
02:30					250				550		
02:45					275	275					
03:00			150		300	300					
03:30						360					
04:00						400					
04:30						450					
05:00			250			500					
05:30			275			550					