

Bombas de infusión SMARTeZ

Infusiones de larga duración



Tabla del efecto de una mayor o menor cantidad de relleno

Nº de modelo (REF)	480050	480080	480010	480030	480250	480090	480260	481200	481170	481180
Caudal (mL/h)	2	2	5	5	10	5	10	20	20	40
Volumen de llenado / reposición (mL)	60	100	60	125	120	270	240	250	500	500
Volumen mínimo (ml)	50	75	60	90	90	250	200	200	350	240
Volumen máximo (ml)	60	120	100	175	170	360	300	300	500	500
Volumen residual (ml)	2	2	2	2	480250	3			5	5

Dauer der Infusion

Horas	Días	Volumen de llenado/reposición (mL)								
6 H										240
8 H						80				320
10 H						100		200		400
12 H			60		120			250	250	500
15 H								300		
18 H			90	90					360	
20 H			100				200		420	
25 H	1 Día	50		125			240		500	
27 H										
30 H		60					300			
36 H	1,5 Días		75	175						
50 H	2 Días		100			250				
54 H			108			270				
72 H	3 Días					360				